


Product Datasheet

Rat Iapp ELISA Kit (orb566889)

Description	Rat Iapp ELISA Kit
Reactivity	Rat
Range	78-5000pg/ml
Target	Iapp
Storage	2-8 °C for 6 months
Note	For research use only
Application notes	Competitive ELISA, Coated with Antigen
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	P12969
Sensitivity	46.9pg/ml
Expiration Date	Please enquire.



Rat Iapp(Islet amyloid polypeptide) ELIS...