


Product Datasheet

Rat ACTH ELISA Kit (orb566894)

Description	Rat ACTH ELISA Kit
Reactivity	Rat
Range	15.625-1000pg/ml
Target	ACTH
Storage	2-8 °C for 6 months
Note	For research use only
Application notes	Competitive ELISA, Coated with Antigen
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Rat
ACTH(Adrenocorticotrophic
Hormone) EL...