


Product Datasheet

Rat Trh ELISA Kit (orb566895)

Description	Rat Trh ELISA Kit
Reactivity	Rat
Range	0.156-10uIU/ml
Target	TRH
Storage	2 - 8 °C for 12 months
Note	For research use only
Application notes	Sandwich ELISA, Double Antibody
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	P01150
Sensitivity	0.094uIU/ml
Expiration Date	Please enquire.



Rat Trh(Prothyroliberin) ELISA Kit

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)