


Product Datasheet

Rat Trh ELISA Kit (orb566895)

| | |
|--------------------------|--|
| Description | Rat Trh ELISA Kit |
| Reactivity | Rat |
| Range | 0.156-10uIU/ml |
| Target | TRH |
| Storage | 2-8 °C for 6 months |
| Note | For research use only |
| Application notes | Sandwich ELISA, Double Antibody |
| Sample Types | Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples |
| Uniprot ID | P01150 |
| Sensitivity | 0.094uIU/ml |
| Expiration Date | Please enquire. |



Rat
Trh(Prothyroliberin)
ELISA Kit