


Product Datasheet

Rat PTH ELISA Kit (orb566898)

Description	Rat PTH ELISA Kit
Reactivity	Rat
Range	12.5-800pg/ml
Target	Pth
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	P04089
Sensitivity	7.5pg/ml
Expiration Date	Please enquire.



Rat
Pth(Parathyroid
hormone)
ELISA Kit