


Product Datasheet

Rat GLP-1 ELISA Kit (orb566923)

Description	Rat GLP-1 ELISA Kit
Reactivity	Rat
Range	0.156-10ng/ml
Target	Glp1r
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	P32301
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Rat
Glp1r(Glucagon-
like peptide 1
recept...