


Product Datasheet

Rat CYP27B1 ELISA Kit (orb566925)

Description	Rat CYP27B1 ELISA Kit
Reactivity	Rat
Range	0.625-40ng/ml
Target	Cyp27b1
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	035132
Sensitivity	0.375ng/ml
Expiration Date	Please enquire.



Rat
Cyp27b1(25-
hydroxyvitamin
D-1 alpha ...