


Product Datasheet

Rat MT-1 ELISA Kit (orb566947)

Description	Rat MT-1 ELISA Kit
Reactivity	Rat
Range	0.312-20ng/ml
Target	MT1
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	P02803
Sensitivity	0.188ng/ml
Expiration Date	Please enquire.



Rat MT-1 (Metallothionein-1) ELISA Kit