


Product Datasheet

Rat Ppargc1a ELISA Kit (orb566955)

Description	Rat Ppargc1a ELISA Kit
Reactivity	Rat
Range	0.156-10ng/ml
Target	Ppargc1a
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	Q9QYK2
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Rat
Ppargc1a(Peroxisome
proliferator-act...