


Product Datasheet

Rat Nts ELISA Kit (orb566965)

Description	Rat Nts ELISA Kit
Reactivity	Rat
Range	1.56-100pg/ml
Target	Nts
Storage	2-8 °C for 6 months
Note	For research use only
Application notes	Competitive ELISA, Coated with Antigen
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	P20068
Sensitivity	0.938pg/ml
Expiration Date	Please enquire.



Rat
Nts(Neurotensin/neuromedin
N) ELISA ...