


Product Datasheet

Rat Chga ELISA Kit (orb566967)

Description	Rat Chga ELISA Kit
Reactivity	Rat
Range	31.25-2000pg/ml
Target	Chga
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	P10354
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Rat
Chga(Chromogranin-
A) ELISA Kit