


Product Datasheet

Rat Gc ELISA Kit (orb566986)

Description	Rat Gc ELISA Kit
Reactivity	Rat
Range	62.5-4000pg/ml
Target	Gc
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	P06883
Sensitivity	37.5pg/ml
Expiration Date	Please enquire.



Rat
Gc(Glucagon)
ELISA Kit