


Product Datasheet

Rat Pgr ELISA Kit (orb566990)

Description	Rat Pgr ELISA Kit
Reactivity	Rat
Range	0.313-20ng/ml
Target	PGR
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	Q63449
Sensitivity	0.188ng/ml
Expiration Date	Please enquire.



Rat
Pgr(Progesterone
receptor) ELISA
Kit