


Product Datasheet

Rat Drd1 ELISA Kit (orb566996)

Description	Rat Drd1 ELISA Kit
Reactivity	Rat
Range	12.5-800ng/ml
Target	Drd1
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	P18901
Sensitivity	7.5ng/ml
Expiration Date	Please enquire.



Rat
Drd1(D(1A)
dopamine
receptor)
ELISA ...