


Product Datasheet

Rat Ggt1 ELISA Kit (orb567012)

Description	Rat Ggt1 ELISA Kit
Reactivity	Rat
Range	1.56-100ng/ml
Target	Ggt1
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	P07314
Sensitivity	0.938ng/ml
Expiration Date	Please enquire.



Rat Ggt1(Gamma-glutamyltranspeptidase 1)...