


## Product Datasheet

### Rat Lbp ELISA Kit (orb567017)

<b>Description</b>	Rat Lbp ELISA Kit
<b>Reactivity</b>	Rat
<b>Range</b>	62.5-4000pg/ml
<b>Target</b>	Lbp
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Sample Types</b>	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
<b>Uniprot ID</b>	<b>Q63313</b>
<b>Sensitivity</b>	37.5pg/ml
<b>Expiration Date</b>	Please enquire.



Rat  
Lbp(Lipopolysaccharide-  
binding prote...