


Product Datasheet

Rat Cyp3a2 ELISA Kit (orb567041)

Description	Rat Cyp3a2 ELISA Kit
Reactivity	Rat
Range	0.625-40ng/ml
Target	Cyp3a2
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	P05183
Sensitivity	0.375ng/ml
Expiration Date	Please enquire.



Rat
Cyp3a2(Cytochrome
P450 3A2) ELISA
Ki...