


Product Datasheet

Rat Cdk5r1 ELISA Kit (orb567047)

Description	Rat Cdk5r1 ELISA Kit
Reactivity	Rat
Range	0.156-10ng/ml
Target	Cdk5r1
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	P61810
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Rat
Cdk5r1(Cyclin-
dependent
kinase 5 act...