


## Product Datasheet

### Rat S100A8 ELISA Kit (orb567075)

<b>Description</b>	Rat S100A8 ELISA Kit
<b>Reactivity</b>	Rat
<b>Range</b>	0.313-20ng/ml
<b>Target</b>	S100A8
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Sample Types</b>	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
<b>Uniprot ID</b>	<b>P50115</b>
<b>Sensitivity</b>	0.188ng/ml
<b>Expiration Date</b>	Please enquire.



Rat  
S100A8(Protein  
S100-A8) ELISA  
Kit