


Product Datasheet

Rat Calprotectin ELISA Kit (orb567076)

Description	Rat Calprotectin ELISA Kit
Reactivity	Rat
Range	0.156-10ng/ml
Target	S100A9
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	P50116
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Rat
S100A9(Protein
S100-A9) ELISA
Kit