


Product Datasheet

Rat Reg4 ELISA Kit (orb567079)

Description	Rat Reg4 ELISA Kit
Reactivity	Rat
Range	1.56-100ng/ml
Target	Reg4
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	Q68AX7
Sensitivity	0.938ng/ml
Expiration Date	Please enquire.



Rat
Reg4(Regenerating
islet-derived prot...