


Product Datasheet

Rat PDGFR beta ELISA Kit (orb567103)

| | |
|------------------------|--|
| Description | Rat PDGFR beta ELISA Kit |
| Reactivity | Rat |
| Range | 0.312-20ng/ml |
| Target | PDGFRb |
| Storage | 2-8 °C for 6 months |
| Note | For research use only |
| Sample Types | Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples |
| Uniprot ID | Q05030 |
| Sensitivity | 0.188ng/ml |
| Expiration Date | Please enquire. |



Rat PDGFRb(Beta-type platelet-derived gr...