


Product Datasheet

Rat Ghrl ELISA Kit (orb567114)

Description	Rat Ghrl ELISA Kit
Reactivity	Rat
Range	4.688-300pg/ml
Target	Ghrl
Storage	2-8 °C for 6 months
Note	For research use only
Application notes	Competitive ELISA, Coated with Antigen
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	Q9QYH7
Sensitivity	2.813pg/ml
Expiration Date	Please enquire.



Rat
Ghrl(Appetite-
regulating
hormone)
EL...