


Product Datasheet

Rat Tfr2 ELISA Kit (orb567129)

Description	Rat Tfr2 ELISA Kit
Reactivity	Rat
Range	0.313-20ng/ml
Target	Tfr2
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	B2GUY2
Sensitivity	0.188ng/ml
Expiration Date	Please enquire.



Rat
Tfr2(Transferrin
receptor
protein 2)...