


Product Datasheet

Rat Nps ELISA Kit (orb567149)

Description	Rat Nps ELISA Kit
Reactivity	Rat
Range	0.156-10ng/ml
Target	Nps
Storage	2-8 °C for 6 months
Note	For research use only
Application notes	Competitive ELISA, Coated with Antigen
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	POC0P7
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Rat
Nps(Neuropeptide
S) ELISA Kit