


Product Datasheet

Rat MDR1 ELISA Kit (orb567154)

Description	Rat MDR1 ELISA Kit
Reactivity	Rat
Range	0.312-20ng/ml
Target	Abcb1
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	P43245
Sensitivity	0.188ng/ml
Expiration Date	Please enquire.



Rat
Abcb1(Multidrug
resistance
protein 1...