


Product Datasheet

Rat ANP ELISA Kit (orb567204)

Description	Rat ANP ELISA Kit
Reactivity	Rat
Range	15.625-1000pg/ml
Target	ANP
Storage	2-8 °C for 6 months
Note	For research use only
Application notes	Competitive ELISA, Coated with Antigen
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	P01161
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Rat
ANP(Atrial
Natriuretic
Peptide)
ELIS...