


Product Datasheet

Rat AQP-2 ELISA Kit (orb567208)

Description	Rat AQP-2 ELISA Kit
Reactivity	Rat
Range	1.563-100ng/ml
Target	AQP-2
Storage	2-8 °C for 6 months
Note	For research use only
Application notes	Competitive ELISA, Coated with Antigen
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	P34080
Sensitivity	0.938ng/ml
Expiration Date	Please enquire.



Rat AQP-2 (Aquaporin 2) ELISA Kit