


Product Datasheet

Rat BMPR-1A ELISA Kit (orb567236)

| | |
|------------------------|--|
| Description | Rat BMPR-1A ELISA Kit |
| Reactivity | Rat |
| Range | 0.781-50ng/ml |
| Target | BMPR-1A |
| Storage | 2-8 °C for 6 months |
| Note | For research use only |
| Sample Types | Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples |
| Uniprot ID | Q78EA7 |
| Sensitivity | 0.469ng/ml |
| Expiration Date | Please enquire. |



Rat BMPR-1A(Bone Morphogenetic Protein R...