


## Product Datasheet

### Rat C1INH ELISA Kit (orb567241)

|                        |  |
|------------------------|--|
| <b>Description</b>     | Rat C1INH ELISA Kit  |
| <b>Reactivity</b>      | Rat  |
| <b>Range</b>           | 7.813-500ng/ml   |
| <b>Target</b>          | C1INH  |
| <b>Storage</b>         | 2-8 °C for 6 months  |
| <b>Note</b>            | For research use only  |
| <b>Sample Types</b>    | Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples |
| <b>Uniprot ID</b>      | <b>Q6P734</b>  |
| <b>Sensitivity</b>     | 4.688ng/ml   |
| <b>Expiration Date</b> | Please enquire.  |



Rat  
C1INH(Complement  
1 Inhibitor) ELISA ...