


Product Datasheet

Rat CASP14 ELISA Kit (orb567257)

Description	Rat CASP14 ELISA Kit
Reactivity	Rat
Range	0.313-20ng/ml
Target	CASP14
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	D3ZZ65
Sensitivity	0.188ng/ml
Expiration Date	Please enquire.



Rat
CASP14(Caspase
14) ELISA Kit