


Product Datasheet

Rat CASP7 ELISA Kit (orb567258)

Description	Rat CASP7 ELISA Kit
Reactivity	Rat
Range	0.313-20ng/ml
Target	CASP7
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	O88550
Sensitivity	0.188ng/ml
Expiration Date	Please enquire.



Rat
CASP7(Caspase
7) ELISA Kit