


Product Datasheet

Rat CDK5 ELISA Kit (orb567275)

Description	Rat CDK5 ELISA Kit
Reactivity	Rat
Range	0.313-20ng/ml
Target	CDK5
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	Q03114
Sensitivity	0.188ng/ml
Expiration Date	Please enquire.



Rat
CDK5(Cyclin
Dependent
Kinase 5)
ELIS...