


Product Datasheet

Rat CFB ELISA Kit (orb567278)

Description	Rat CFB ELISA Kit
Reactivity	Rat
Range	3.125-200ng/ml
Target	CFB
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	Q6MG74
Sensitivity	1.875ng/ml
Expiration Date	Please enquire.



Rat
CFB(Complement
Factor B) ELISA
Kit