


Product Datasheet

Rat CORT ELISA Kit (orb567302)

Description	Rat CORT ELISA Kit
Reactivity	Rat
Range	2.81-180ng/ml
Target	CORT
Storage	2-8 °C for 6 months
Note	For research use only
Application notes	Competitive ELISA, Coated with Antigen
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Sensitivity	1.688ng/ml
Expiration Date	Please enquire.



Rat
CORT(Corticosterone)
ELISA Kit