


Product Datasheet

Rat CP ELISA Kit (orb567303)

Description	Rat CP ELISA Kit
Reactivity	Rat
Range	6.25-400ng/ml
Target	CP
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	P13635
Sensitivity	3.75ng/ml
Expiration Date	Please enquire.



Rat
CP(Ceruloplasmin)
ELISA Kit