


Product Datasheet

Rat CTXI ELISA Kit (orb567317)

Description	Rat CTXI ELISA Kit
Reactivity	Rat
Range	0.156-10ng/ml
Target	CTXI
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Rat
CTXI(Cross
Linked C-
telopeptide
of T...