


Product Datasheet

Rat DBP ELISA Kit (orb567333)

Description	Rat DBP ELISA Kit
Reactivity	Rat
Range	15.625-1000ng/ml
Target	DBP
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	P04276
Sensitivity	9.375ng/ml
Expiration Date	Please enquire.



Rat
DBP(Vitamin
D Binding
Protein)
ELISA...