


Product Datasheet

Rat ENA-78 ELISA Kit (orb567353)

Description	Rat ENA-78 ELISA Kit
Reactivity	Rat
Range	15.625-1000pg/ml
Target	ENA-78
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	P97885
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Rat ENA-78(Epithelial Neutrophil Activat...