


Product Datasheet

Rat EP2 ELISA Kit (orb567354)

Description	Rat EP2 ELISA Kit
Reactivity	Rat
Range	39.063-2500pg/ml
Target	EP2
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	Q62928
Sensitivity	23.438pg/ml
Expiration Date	Please enquire.



Rat
EP2(Prostaglandin
E Receptor 2)
ELIS...