


Product Datasheet

Rat FASL ELISA Kit (orb567370)

Description	Rat FASL ELISA Kit
Reactivity	Rat
Range	31.25-2000pg/ml
Target	FASL
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	P36940
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Rat
FASL(Factor
Related
Apoptosis
Ligand...