


Product Datasheet

Rat FASN ELISA Kit (orb567371)

Description	Rat FASN ELISA Kit
Reactivity	Rat
Range	0.313-20ng/ml
Target	FASN
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	P12785
Sensitivity	0.188ng/ml
Expiration Date	Please enquire.



Rat
FASN(Fatty
Acid
Synthase)
ELISA Kit