


Product Datasheet

Rat FGF-10 ELISA Kit (orb567374)

| | |
|------------------------|--|
| Description | Rat FGF-10 ELISA Kit |
| Reactivity | Rat |
| Range | 15.625-1000pg/ml |
| Target | FGF-10 |
| Storage | 2-8 °C for 6 months |
| Note | For research use only |
| Sample Types | Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples |
| Uniprot ID | P70492 |
| Sensitivity | 9.375pg/ml |
| Expiration Date | Please enquire. |



Rat FGF-10(Fibroblast Growth Factor-10) ...