


Product Datasheet

Rat FGF21 ELISA Kit (orb567376)

Description	Rat FGF21 ELISA Kit
Reactivity	Rat
Range	31.25-2000pg/ml
Target	FGF21
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	Q8VI80
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Rat
FGF21(Fibroblast
Growth Factor
21) E...