


Product Datasheet

Rat FGFR1 ELISA Kit (orb567379)

Description	Rat FGFR1 ELISA Kit
Reactivity	Rat
Range	0.156-10ng/ml
Target	FGFR1
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	Q04589
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Rat
FGFR1(Fibroblast
Growth Factor
Recep...