


Product Datasheet

Rat FSH ELISA Kit (orb567384)

Description	Rat FSH ELISA Kit
Reactivity	Rat
Range	2.344-150mIU/ml
Target	FSH
Storage	2-8 °C for 6 months
Note	For research use only
Application notes	Competitive ELISA, Coated with Antigen
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	P18427
Sensitivity	1.406mIU/ml
Expiration Date	Please enquire.



Rat
FSH(Follicle
Stimulating
Hormone)
EL...