


Product Datasheet

Rat GAL1 ELISA Kit (orb567393)

Description	Rat GAL1 ELISA Kit
Reactivity	Rat
Range	0.625-40ng/ml
Target	GAL1
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	P11762
Sensitivity	0.375ng/ml
Expiration Date	Please enquire.



Rat
GAL1(Galectin
1) ELISA Kit