


Product Datasheet

Rat GCP-2 ELISA Kit (orb567404)

| | |
|------------------------|--|
| Description | Rat GCP-2 ELISA Kit |
| Reactivity | Rat |
| Range | 31.25-2000pg/ml |
| Target | GCP-2 |
| Storage | 2-8 °C for 6 months |
| Note | For research use only |
| Sample Types | Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples |
| Uniprot ID | G3V6C8 |
| Sensitivity | 18.75pg/ml |
| Expiration Date | Please enquire. |



Rat GCP-
2(Granulocyte
Chemotactic
Protei...