


Product Datasheet

Rat GDF15 ELISA Kit (orb567408)

Description	Rat GDF15 ELISA Kit
Reactivity	Rat
Range	31.25-2000pg/ml
Target	GDF15
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	Q9Z0J6
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Rat
GDF15(Growth
Differentiation
Factor ...