


Product Datasheet

Rat GH ELISA Kit (orb567416)

Description	Rat GH ELISA Kit
Reactivity	Rat
Range	9.375-600pg/ml
Target	GH
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, Cell or tissue lysate, Other liquid samples
Uniprot ID	P01244
Sensitivity	5.625pg/ml
Expiration Date	Please enquire.



Rat
GH(Growth
Hormone)
ELISA Kit